

# GMI ELECTROSURGICAL PEN DR/DB - SURGICAL LINE




## FREEDOM

- Universal three-pin connector.
- 10-foot cable.


## SAFETY

- Actuation Buttons
- *Push Button* for easy handling.
- Effective hemostasis.
- Precise cut.



## VARIETY

- Also available in the *Rocker Switch* version. (sale on request)



THE GMI  
ELECTROSURGICAL  
PEN, WAS  
DESIGNED TO FIT  
PERFECTLY IN THE  
HANDS, AVOIDING  
FATIGUE.

**COMFORT**  
- Anatomical design to  
reduce hand strain during  
use.



### INDICATION

The electro-surgical pen is a versatile and efficient device. The precise activation facilitates its use during the surgical procedure. When in operation, the pen uses an alternating current, raising the temperature of the tissue to perform functions such as: cutting, coagulation and dissection.

DISPOSABLE ELECTROSURGICAL PEN GMI ANVISA Registration No. 80423540045				
Code	ACTUATION TYPE	BLADE LENGTH (cm)	BLADE MATERIAL	Cable Length (cm)
DB	Push button	2,5cm	Stainless Steel	300
DR (Sale on request)	Rocker Switch	2,5cm	Stainless Steel	300