


### VERSATILITY

- Gastrostomy probe with traction loop in several sizes.



### RELIABILITY

- The items in the kit meet the needs of passing and fixing the gastrostomy probe, in a safe and effective way

THE PULL  
GASTROSTOMY  
KIT IS A  
COMPLETE,  
QUALITY OPTION  
FOR THE  
PRACTICE MOST  
USED BY  
PROFESSIONALS.



COMPLETE KIT TO PERFORM THE  
PROCEDURE

INDICATION

The kit for percutaneous endoscopic gastrostomy by the PULL method, aims to obtain an access route for prolonged enteral nutrition, in cases of patients with acute or chronic clinical conditions that present malnutrition or nutritional risk.

INCLUDED ITEMS

- 01 Gastrostomy Probe with traction handle.
- 01 Puncture needle.
- 01 Polypectomy forceps.
- 01 Fixing ring.
- 01 Double connector.
- 01 Enfit double connector.
- 01 Surgical field.
- 01 10ml syringe.
- 01 No. 11 Scalpel
- 01 Flexible guide.
- 01 Clamp

	355   GMI PEG PULL ENDOSCOPIC GASTROSTOMY KIT			
	ANVISA Registration No. 80423540034			
Code	FR	PROBE DIAMETER (mm)	NEEDLE DIAMETER (G)	Probe Length (cm)
355-16	16	5,3	16	60
355-18	18	6,0	16	60
355-20	20	6,6	14	60
355-24	24	8,0	14	60