

# GMI PERCUTANEOUS GASTROSTOMY BUTTON

## 350 - PERCUTANEOUS GASTRIC DIGESTIVE LINE



### COMFORT

- Three diameter options.
- Anti-reflux valve.
- Medical-grade silicone.



### SAFETY

- Safe and stable fixation.
- High durability and maximum comfort.



### PRACTICALITY

- Kit with indispensable components for the application of the button, including tube for feeding and gastric decompression.

THE GMI  
GASTROSTOMY  
BUTTON WAS  
DEVELOPED  
TO ENABLE  
GREATER  
COMFORT TO  
THE PATIENT IN  
THE  
GASTROSTOMY.  
HIGHLY  
DURABLE, THE  
DEVICE  
PROMOTES  
SAFETY AND  
WELL-BEING.



#### KIT COMPONENTS

- 01 Button insertion device.
- 01 Abdominal wall meter.
- 01 60ml syringe.
- 01 5ml syringe.

#### MULTIPLE CONNECTIONS

- 02 Y-extender tubes with clamp.
- 01 Extender tube for gastric decompression.

The accessory may be marketed separately:

GMI TUBE	
ANVISA Registration No. 80423540016	
CODE	FR
TE-350-03*1	12
TE-350-03*2	17
TE-350-03*3	24
TE-350-03*4	34
TE-350-03*5	44

350   GMI PERCUTANEOUS GASTROSTOMY BUTTON				
ANVISA Registration No. 80423540024				
Code	FR	EXTERNAL DIAMETER (mm)	INTERNAL DIAMETER (mm)	Length (cm)
350-18-12	18	6.0	3.5	1.2
350-18-17	18	6.0	3.5	1.7
350-18-24	18	6.0	3.5	2.4
350-18-34	18	6.0	3.5	3.4
350-20-17	20	6.7	4.0	1.7
350-20-24	20	6.7	4.0	2.4
350-24-12	24	8.0	5.0	1.2
350-24-17	24	8.0	5.0	1.7
350-24-24	24	8.0	5.0	2.4
350-24-34	24	8.0	5.0	3.4
350-24-44	24	8.0	5.0	4.4