

PRACTICALITY

- Y-connector that facilitates handling.
- Radiopaque to the fullest extent.

BIOCOMPATIBLE


- Developed in medical grade silicone.

COMFORT

- Atraumatic finish that facilitates passage.

DURABILITY

- Silicone connector resistant to deformation and rupture.



THE GMI
NASOGASTRIC
C TUBE,
ALLOWS EASY
INSERTION
AND PATIENT
COMFORT
WITH THE
ADVANTAGES
OF MEDICAL
GRADE
SILICONE.

VERSATILITY

- Conexión en Y que facilita el montaje del equipo



INDICATION

- The GMI Nasogastric Probe was developed to meet all patient profiles, especially those requiring emptying or gastric drainage.

CHARACTERISTICS

- Y-connection with lid.
- Graduated every 05cm and marked every 01cm.
- Radiopaque line.
- Sterile.

	391/392 GMI NASOGASTRIC PROBE ANVISA Registration No. 80423540005			
Code	FR	EXTERNAL DIAMETER (mm)	INTERNAL DIAMETER (mm)	Length (cm)
391-10	10	3.30	2.00	105
391-12	12	4.00	2.50	120
392-04	04	1.35	0.75	45
392-06	06	2.00	1.00	45
392-07	07	2.30	1.25	45
392-09	09	3.00	1.75	45