

**COMPLETE KIT FOR
VENTILATORY SUPPORT**

Extremely versatile, with adapters that allow use in various mechanical fans available on the market.

ACCURACY

Pressure monitoring line for complete and reliable data.



COMFORT

Extra soft cotton cap, ensuring comfort for the newborn and with adjustable straps for perfect tube fixation.

ANATOMICAL

100% silicone nasal prong, with delicate adjustment to the neonate's skin.

**GMI COMPLETE
NEONATAL
NASAL C.P.A.P.
SYSTEM, IN
100% SILICONE,
GUARANTEES
AN EXCELLENT
ANATOMICAL
ADAPTATION
EVEN IN
EXTREMELY
LOW WEIGHT
NEONATES.**



VERSATILITY

Complete kit with adapters that allows the use in several brands of mechanical ventilators.

INDICATION

Use specially designed to attend neonatology and pediatrics, to provide continuous positive air pressure, for ventilatory support and flexibility necessary for daily clinical routines.

The GMI C.P.A.P. allows the use in several models of mechanical ventilators, due to the versatility of the adapters included in the kit.

ITEMS INCLUDED

- Complete system for neonatal nasal CPAP:
- 01 Nasal prong in 100% silicone with anatomical fit.
 - 02 120cm corrugated tubes for respiratory branch. 01 40cm corrugated tube for humidifier.
 - 01 Pressure monitoring line.
 - 01 Extra soft cotton cap with adjustable velcros for tube stabilization.
 - 01 Adapter for silicone respirator.
 - 02 90° connectors for inspiratory/spiratorial branch.
 - 02 22mm Connectors.
 - 03 10mm Standard Connectors.

233 GMI NASAL C.P.A.P. ANVISA Registration No. 80423540003					
Code	CAP SIZE	FR	CANNULA DIAMETER (mm)	CANNULA LENGTH (mm)	Recommended Weight (g)
233-00	P	08	2,90	8,5	Up to 1000
233-01	P	10	3,30	8,8	1000-1500
233-02	M	12	4,15	10,0	1500-2000
233-03	G	15	4,92	12,0	more than 2000